

Development of a custom CMOS modulated light camera and its application to picosecond ultrasonics

Richard Smith, Roger Light,
Mike Somekh, Steve Sharples,
Mark Pitter, Nick Johnston

Applied Optics Group

Electrical Systems and Optics Research
Division

University of Nottingham





The University of
Nottingham

Talk Outline

- Introduction
- Custom detector
- Coating thickness
- 2D thickness maps
- Future work
- Conclusions



Introduction – Laser generation and detection of ultrasound

- Laser pulse absorbed
- Rapid local heating
- Heating causes expansion
- Expansion generates sound wave
- Different mechanisms for detection
- Reflectivity
- Surface changes
- Very large background with small signal of interest



Experiment Setup

- Typical pump/probe setup
- Time delay between pump and probe imposed by mechanical scan of delay line mirror
- Pump beam modulated by mechanical chopper
- Pump and probe beams separated by polarisation optics
- Detection usually consists of photodiode and lock in amplifier



Multiple channel detection

- Multiple lock-ins become impractical – use phase stepping to demodulate signal
- Capture enough photons to achieve required SNR – use integrating detectors with large well depth
- Be aware of amplitude phase cross talk issues when using squarewave modulation

$$Phase = \tan^{-1} \left(\frac{S_2}{S_1} \right)$$

$$Amplitude = \sqrt{S_1^2 + S_2^2}$$

α = angle

$\Delta\alpha$ = change in α between steps

S_1 = real part of signal

S_2 = imaginary part of signal

N = number of steps


m = current step

I = measured intensity



$$\Delta\alpha = \frac{2\pi}{N}$$

$$S_1 = \sum_{m=0}^{N-1} I_m \cos \alpha_m$$

$$S_2 = \sum_{m=0}^{N-1} I_m \sin \alpha_m$$





Custom CMOS detector



- 64x1 linear array detector
- Pixels built on active sensor principle with 4 large independently switchable capacitors to increase well depth
- 4 phase mode of operation : reset, integration, idle and readout
- Global shutter removes the phase shift between pixels caused by the rolling shutter in commercial detector
- Pixels are randomly addressable
- Faster readout (frame rate of 160KHz /10MHz pixel rate)



Spatial information

- Cylindrical lenses added to produce line focus for both pump and probe beams
- Each pixel now corresponds to a different spatial location on the surface of the sample
- Increase the overall optical power used in the system as the light is spread out over an extended region





Layer Thickness

- Pump light is absorbed by the top layer
- Acoustic pulse travels through the layer and is partially reflected by interface
- When pulse arrives back at top surface it changes the amount of probe light reflected back to the detector




Layer Thickness

- Chromium on silicon sample
- First 3 echoes are visible
- ~30 times faster than single point detector in this case (112 seconds for array data)
- Signals at extremities are smaller due to reduced light intensities
- Sample has two regions of different thickness
- Echoes from the different regions and the transition between them clearly visible
- Echo locations can be converted into thicknesses if the velocity is known




2D Thickness maps

- Using the first echo location we can use fewer averages
- This allows the possibility of mapping a large area of the sample in a practical amount of time.
- Top figure took a couple of hours and shows a region of the sample where the thicker layer has a region of damage (probably a surface scratch)
- The bottom figure shows a composite image mapping regions either side and across the transition
- The results on the flat region are still good even with the reduced averaging– the standard deviation of the thickness in the inset is $\sim 0.25\text{nm}$
- This image took $\sim 2.45\text{hrs}$ to acquire




Future Work

- The current limit to speed of this system is the ADC and delay stage we employ.
- An additional speed increase of at least $\sim 5x$ is possible if these are replaced.
- Improve base SNR by looking into laser noise cancellation
- Frequency double one beam to allow more complex detection schemes and completely remove pump laser noise.



Conclusions

- Our detector is ~ 30 times faster than our single point detector system
- We can now perform measurements such as thickness maps on thin coatings over a large area of sample
- Additional speed increases of ~ 5 times should be easy to achieve



- Any Questions?

Acknowledgements:

We wish to thank the EPSRC for funding this work